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General Regulations of the laboratories of experimental psychology (Experimental Psychology Lab) , CIMeC - Fedrigotti Palace, Rovereto

December 2013

Rules and Procedures for the use of the EEG laboratory

Access to the laboratory

- The coordinator of the laboratories shall establish, through advice from expert users (independent users), those users who are autonomous in the use of the laboratory, and make sure that all users are aware of the training procedures and standards for use of the laboratory.
- It is important to plan well in advance to use the EEG laboratory. For example, if it is necessary that the experimenter is trained about the procedures for mounting, washing and disinfecting the EEG caps, the experimenter is required to attend the laboratory EEG (assisting other independent users during the experimental sessions) a few months in advance before the starting of the experiment.
- The independent users identified by the lab coordinator, will test whether an experimenter can access the independent user-level by assisting an experimental session in which the experimenter proves to be autonomous in handling the phases of affixing the electrodes, recording and cleaning the caps.

General rules of use

- All equipment (caps, electrodes, preparation room and recording booth) must be kept clean and tidy.
- You cannot change the standard electrode montages (see Annex A) without permission of the laboratory technician. For efficient management of procedures for affixing the electrodes, it is also forbidden to move the consumables (stored in the closets of the preparation and recording rooms).
- Users are required to notify the lab technician when the consumable material is about to end. Information on the amount of consumables required for each experimental session are listed in annex B.
- If a user suspects a malfunction of the hardware and software of the laboratory (for example, difficulty in measuring impedance and / or EEG recording , the presence of unusual electrical artifacts), he/she is required to promptly notify the lab technician and to take note of these malfunctions (e.g, by reporting which cap, configuration, and connectors you were using at the time of malfunction).

CiMeC

Palazzo Fedrigotti
corso Bettini, 31 - 38068 Rovereto, Trento (Italy)
ph. +39 0464 808615, fax +39 0464 808690
e-mail: cimec@unitn.it
<http://www.cimec.unitn.it>

- Each user is required to create a separate workspace for recording EEG data, and to make sure not to change the workspace of other users. In addition, any changes in the " digital port setting" should be reset to standard after use.

Procedures for affixing, washing and disinfection of electrodes

- Mounting "EasyCap" electrode caps. For a correct positioning of the caps, you should refer to the information contained in the manual EasyCap, located in the preparation room and available in English and Italian. All procedures listed in the manual for mounting the caps should be followed in detail. During all phases of affixing the electrodes and caps users should wear disposable latex gloves and wash their hands before preparing the consumables
- • Connecting the electrodes to the amplifiers and the recording system. The cap adapter/s needs to be connected to the adapter located in the main input-box (the element connected to the amplifiers) with the utmost care, as the input box is very delicate and easily wearable. Care must also be observed in the case of direct connection of the electrodes (eg, ground and reference) in the inputbox. The ground electrode must be connected to the input-box first.
- Cleaning and disinfection of the caps. All procedures listed in this phase are performed wearing gloves. After removing the caps from the subject, you should immediately wash and disinfect the caps and the electrodes. If this is not possible, because for example the participant want to wash their hair, the cap should be soaked in warm water and mild shampoo until you can begin the proper washing. Make sure that the adapters and the pins of the electrodes do not get in contact with water.
- First, wash the caps (internally and externally) with warm water jet and/or with the electric toothbrush, so as to eliminate most of the remaining gel from the electrodes. Second , brush all electrodes (internally and externally) with the toothbrush and the shampoo or soap, until there is no longer any trace of gel.
Once rinsed with warm water, the cap must be put in disinfectant (Metricide) and left in the disinfecting solution for at least 30 minutes. Pay attention to the expiration date of the disinfectant (approximately 4 weeks after activation). If the disinfectant is expired, you must notify the technician. After the time for disinfection, thoroughly rinse the cap with warm water and place it to dry on the box on the washing machine , always making sure that the plugs of the electrodes and the adapter never get in contact with water or wet material.
- The caps should never be left to soak for longer times than those required for cleaning and disinfection.
- All the consumables used (gloves, disposable syringes, cotton swabs, plastic spoons, paper tissues, plastic stickers, etc) must be eliminated at the end of the experimental session before a new participant enters the room. Failure to

follow the procedures listed above may compromise the hygienic conditions of the laboratory, with potential risks for subsequent experimental sessions.

- The towels used during the experiment should be placed in the washing machine and once there are about 8 used towels, a washing and drying cycle can be programmed when there are no running experiments.

Appendix A: Standard montages- EasyCap (EEG laboratory)

As of today two are two standard electrode configurations:

- a configuration with a low density array- CFG1 - with 23 sites (Figure 1) , plus any sites for facial electro-oculogram;
- configuration with higher density -CFG2-, corresponding to the Easycap M11 montage, with 63 sites (Figure 2), plus a maximum of three sites for the electro-oculogram.

If for some studies there is the need to use intermediate or different configurations it is advisable that these non-standard configurations be obtained by changing the settings to a minimum level, as fixing the electrodes into the slots of the caps is time-consuming and risky, because it is possible to inadvertently swap the position of sites.


Figure 1: Mounting CFG1
CFG2


Figure 2: Mounting

There are different caps with size 54cm , 56cm , 58 cm, 60 cm for CFG1 (measures in cm circumference of head 54, 56, 58, 60) and 4 caps CFG2 (sizes 54, 56,58,60).

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If you use different configurations, or an additional number of electrodes for CGF1 , after the last use in a slot reservation you must put the lab material back to the standard configuration; this operation should be done as soon as possible (after the time to wash and disinfect dry material) and no later than one day from the last use. Users using non-standard configurations are required to notify the other users of the laboratory of the exact timing of the reconfiguration of the standard settings.

CiMeC

Palazzo Fedrigotti
corso Bettini, 31 - 38068 Rovereto, Trento (Italy)
ph. +39 0464 808615, fax +39 0464 808690
e-mail: cimec@unitn.it
<http://www.cimec.unitn.it>

Appendix B: Use of consumables per experimental session with EasyCaps

The amount of consumables will vary with the type of montage used. It is assumed that CFG1 is performed by a single operator and CFG2 by two operators.

<i>DESCRIPTION</i>	<i>CFG2</i>	<i>CFG1</i>
<i>Cotton swabs</i>	4	2
<i>Syringes 20cc</i>	2	1
<i>Abrasive gel</i>	40 ml	20 ml
<i>Disposable gloves (pair)</i>	2	1